

Test results of the Nocolyse sporicide activity against *Clostridium difficile* using the ASTM 2197:2002 test method (vapour phase)

Study done by :

Laboratoire du Centre de Référence *Clostridium difficile* belge
Faculté de Médecine de l'Université Catholique de Louvain
Unité de Microbiologie
Avenue Hippocrate, 54 UCL 54.90
B-1200 Bruxelles

Professor M. Delmée : Head of Laboratory
Docteur A. Simon : Hygienist-microbiologist
Mr.Johann Van Broeck : Head of Study

With the collaboration of :

Hôpital Iris-Sud à Bruxelles, Site Bracops à Anderlecht
Rue Dr Huet,79 B-1070 Bruxelles

Dr. J.M.Hubrechts, PhD : Laboratory Director

TABLE OF CONTENTS

TABLE OF CONTENTS.....2

1. Summary of the test3

2. Introduction.....3

 Test Method.....4

3. Test Product.....4

4. Tested Microorganisms.....5

5. Study Facility.....5

6. Study Period.....6

7. Experimental Procedure.....7

 Carriers.....7

 Culture Media and Solutions.....7

 Preparation of the Spores.....7

 Disinfection Testing and Calculating Log₁₀ reductions.....7

8. Results.....9

9. Conclusions.....10

10. Annex 1: Description of the Test Rooms.....12

11. Annex 2: Plan of the Test Rooms.....14

12. Annex 3: Testing Parameters.....16

1. Summary of the test

Tested product: Nocolyse/Nocospray System

Manufacturer: Airel

Method used: Derivative of ASTM E 2197-02

Carriers: Polystyrene petri dishes

Testing temperature: Room temperature ($22 \pm 2^\circ\text{C}$)

Dilution: None required

Four rooms (2 rooms per batch) were studied with 5 positions per room

Results:

Table I: Test results of the Nocolyse sporicide activity against *Clostridium difficile* using the ASTM 2197:2002 test method (vapour phase)

| Organisms Tested | Positions | Room Numbers | | | |
|--------------------------------|-----------|--------------|------|------|------|
| | | 310 | 311 | 304 | 305 |
| <i>C. difficile</i> ATCC 43598 | 1 | 7.00 | 6.00 | 5.37 | 5.63 |
| | 2 | 6.25 | 5.65 | 5.43 | 5.60 |
| | 3 | 6.88 | 5.49 | 5.55 | 5.66 |
| | 4 | 6.78 | 5.43 | 5.72 | 5.65 |
| | 5 | 5.06 | 5.36 | 5.51 | 5.62 |
| <i>C. difficile</i> NAP1A/027 | 1 | 6.52 | 5.00 | 5.60 | 5.55 |
| | 2 | 5.89 | 5.17 | 5.52 | 6.00 |
| | 3 | 6.00 | 6.10 | 6.10 | 5.52 |
| | 4 | 5.74 | 6.10 | 5.74 | 5.64 |
| | 5 | 4.97 | 4.86 | 5.64 | 4.80 |

As the testing was not done in the laboratory, closable petri dishes were chosen in order to facilitate the transportation of the inoculum to and from the hospital test site.

2. Introduction

The objective of the study is to determine the efficacy, in a hospital milieu, of the Nocolyse/Nocospray System against *Clostridium difficile* ATCC 43598 and NAP1A/027.


The Nocolyse/ Nocospray System consists of a portable plug-in machine (Nocospray) and a liquid hard surface disinfectant (Nocolyse). The disinfecting liquid Nocolyse is nebulized at 37°C to form a dry vapour which is then propelled by the Nocospray turbine at 80m/sec permitting a rapid distribution throughout the target space.



The Microbiology testing was conducted at the Microbiology Laboratory of UCL by professor M. Delmee. The field testing was done at the Iris-Sud Hospital, Bracops in Bruxelles with the collaboration of doctor J.M. Hubrechts.

Test Method

The testing of the vapour was done in four standard hospital rooms with petri dishes of each strain located in five different positions. The population of each strain was counted on the carriers and compared to the control carriers not exposed to the Nocolyse/Nocospray System. The test was done using a derivative of the quantitative carrier test method, ASTM E 2197 (2002).



3. Test Product

Name: NOCOLYSE

Manufacturer: Airel, Champigny-Sur-Marne, France

Active Material: 6% solution of hydrogen peroxide.

Principal: Liquid Nocolyse is vaporized at 37°C into a dry vapour within the space to be disinfected.

4. Tested microorganisms

The microbial strains used for the testing are illustrated in Table II

Table II: Microbial strains utilized

| Microorganisms | Strains | Provider | Lot Numbers |
|------------------------------|-----------------------|------------------------|-------------|
| <i>Clostridium difficile</i> | ATCC 43598 (B+/A-) | Prof. M. DELMEE UCL | 1470 |
| <i>Clostridium difficile</i> | NAP1A/027 | M. WARNY ACAMBIS | 1067 |

5. Study Facility

The hospital test site was at the Iris Sud Hospital, Bracops in Bruxelles. Four standard hospital rooms were used for the test (refer to Annexe 1). The carriers were prepared and analyzed at the Microbiology Lab at the Catholique Louvain University.

Strains :

Two strains of *Clostridium difficile* were used, the Nord American epidemic strain NAP1A/027 and strain ATCC 43598 (strain B+/A-).

Preparation of the spores:

Inoculate 10 mL of Brain Heat Infusion broth (BHI) – reference BD BBL 211059 – with a colony of *Clostridium difficile* grown on a Columbia agar 24 hours (BD 017423). After 24 hours of anaerobic incubation, dilute the broth 50 times in BHI and allow the colonies to grow for 10 days anaerobically. Centrifuge to reduce 10 times.

Homogenize well and place 1 mL of the solution in an empty petri dish. Add 9 mL of absolute ethanol, mix well and let it evaporate slowly (laminar flow).

Neutralization:

Even though a neutralization process was not necessary due to the rapid degradation of hydrogen peroxide vapour particles at a temperature of 37°C, a neutralization was done in order to follow lab procedures. A neutralization was done right after exposure of the carriers to the Nocolyse/Nocospray System with 10mL of Lethen broth (polysorbate 80, reference BD 268110).

Spore Count:

The number of spores of the control carrier and the exposed carrier were counted after successive dilutions of 10 in 10 of the neutralized solution in saline water. 200 μL of each dilution were spread on a TCCFA medium (taurocholate-containing medium). After an incubation period of 48 hours at a temperature of $35^{\circ}\text{C} \pm 2^{\circ}\text{C}$, the number of colonies was counted and averaged.

6. Study Period

The test was done on the 28th of July 2009 at the Bracops Hospital.

Five petri dishes containing spores of the NAP1A/027 strain and the ATCC 43598 strain were placed in four different hospital rooms.

The control petri dish ATCC 43598 contained 30 000 000 spores/mL which equals to 1 265 823 spores per cm^2 .

The control petri dish NAP1A/027 contained 10 000 000 spores/mL which equals to 421 940 spores per cm^2 .

In order to accommodate the study, different volumes than the ones stated in the standard were used. The final number of cells was equal to the one stated in the ASTM standard.



7. Experimental Procedure

The testing done was based on a derivative of the ASTM E2197 :2002 standard.

Carriers :

The testing was done on 54mm diameter sterile petri dishes.

Culture Media and solutions :

Brain Heart Infusion broth (BHI) - (BD BBL 211059)

Columbia Agar With 5% Sheep Blood (BD017423)

Absolute ethanol

Lethen broth used with polysorbate 80 (BD 268110)

TCCFA medium (taurocholate-containing medium)

Preparation of the spores :

Inoculate 10 mL of Brain Heat Infusion broth (BHI) – reference BD BBL 211059 – with a colony of *Clostridium difficile* grown on a Columbia agar 24 hours (BD 017423). After 24 hours of anaerobic incubation, dilute the broth 50 times in BHI and allow the colonies to grow for 10 days anaerobically. Centrifuge to reduce 10 times.

Make sure to homogenize well and place 1 mL of the solution in an empty petri dish. Add 9 mL of absolute ethanol, mix well and let it evaporate slowly (laminar flow).

In order to accommodate the study, different volumes than the ones stated in the standard were used. The final number of cells was equal to the one stated in the ASTM standard.

Disinfection Testing

Disinfection testing took place at the Anderlecht Bracops Hospital under the supervision of the doctor J.M. Hubrechts, hygienist doctor of the institution.



5 dishes of each strain were located in various positions in each of the 4 rooms used for the testing. The rooms were standard hospital rooms with a 50m³ volume (refer to Annexe 2).



The petri dish covers were removed prior to disinfection with the Nocolyse/Nocospray System which was placed at a height of 40 cm (as seen in the picture below) giving a wide diffusion angle.



The room volume selector was set at seven times the volume of the room since experience has shown this factor to be effective in killing all microorganisms that have been previously evaluated using 350 mL of Nocolyse. Two lots of Nocolyse were tested (2 rooms per lot). One hour after the Nocospray stopped, neutralization was done using 10 ml of Lethen broth. The petri dishes were closed and sent to the UCL Faculty of Medicine where the analysis was done. The dishes were kept at 4°C during transportation. Three control petri dishes were kept at the Laboratory (non- exposed to the Nocolyse/Nocospray System).

Spore Count :

The number of spores of the control carrier and the exposed carrier were calculated by successive dilutions 10 in 10 in the neutralized solution in saline water. 200 μ L of each dilution were spread on a TCCFA medium (taurocholate-containing medium). After an incubation period of 48 hours at a temperature of 35°C +/- 2°C, the number of colonies was counted and averaged.

Calculating Log₁₀ Reductions

The log₁₀ reduction was calculated for each carrier in relation to the mean of the 3 controls.

Reduction on the Carrier $n = \text{Log}_{10}(T) - \text{Log}_{10}(N_n)$

T= mean of the number of colonies for the controls

N_n= number of colonies on the carrier n (n= 1 to 5). If N_n = 0, we add 1 in order to calculate the log.

The product Nocolyse reached the 6 log kill target and was able to eliminate 99.9999% of the organisms which were tested.

8. Results

Table VII: Log Reduction Results for *Clostridium difficile* ATCC 43598

| Rooms | Lot Numbers | Positions | Spores/mL | Spores/cm ² | T | R. Log |
|-------|-------------|-----------|-----------|------------------------|--------|--------|
| 310 | 060509 OS | 1 | 0 | 0 | 3 E+07 | 7,00 |
| | | 2 | 17 | 1 | 3 E+07 | 6,25 |
| | | 3 | 2 | 0 | 3 E+07 | 6,88 |
| | | 4 | 5 | 0 | 3 E+07 | 6,78 |
| | | 5 | 262 | 11 | 3 E+07 | 5,06 |
| 311 | 220709 OS | 1 | 30 | 1 | 3 E+07 | 6,00 |
| | | 2 | 67 | 3 | 3 E+07 | 5,65 |
| | | 3 | 97 | 4 | 3 E+07 | 5,49 |
| | | 4 | 112 | 5 | 3 E+07 | 5,43 |
| | | 5 | 132 | 6 | 3 E+07 | 5,36 |
| 304 | 060509 OS | 1 | 127 | 5 | 3 E+07 | 5,37 |
| | | 2 | 112 | 5 | 3 E+07 | 5,43 |
| | | 3 | 85 | 4 | 3 E+07 | 5,55 |
| | | 4 | 57 | 2 | 3 E+07 | 5,72 |
| | | 5 | 92 | 4 | 3 E+07 | 5,51 |
| 305 | 220709 OS | 1 | 70 | 3 | 3 E+07 | 5,63 |
| | | 2 | 75 | 3 | 3 E+07 | 5,60 |
| | | 3 | 65 | 3 | 3 E+07 | 5,66 |
| | | 4 | 67 | 3 | 3 E+07 | 5,65 |
| | | 5 | 72 | 3 | 3 E+07 | 5,62 |

T = 5,65 : mean of the spore count on non-exposed carriers (non-exposed to the Nocolyse/Nocospray System)

N = mean of the spore count on exposed carriers placed in different locations in the four rooms

R. Log = Reduction of the strain population by the disinfection system in comparison to the mean of the non exposed carriers (T)

Table VIII: Log Reduction Results for *Clostridium difficile* NAP 1A/027

| Rooms | Lot Numbers | Positions | Spores/mL | Spores/cm ² | T | R. Log |
|-------|-------------|-----------|-----------|------------------------|---|--------|
| 310 | 060509 OS | 1 | 3 | 0 | 1 | 6,52 |
| | | 2 | 13 | 1 | 1 | 5,89 |
| | | 3 | 10 | 0 | 1 | 6,00 |
| | | 4 | 18 | 1 | 1 | 5,74 |
| | | 5 | 108 | 2 | 1 | 4,97 |
| 311 | 220709 OS | 1 | 100 | 4 | 1 | 5,00 |
| | | 2 | 68 | 3 | 1 | 5,17 |
| | | 3 | 8 | 0 | 1 | 6,10 |
| | | 4 | 8 | 0 | 1 | 6,10 |
| | | 5 | 138 | 6 | 1 | 4,86 |
| | | 1 | 25 | 1 | 1 | 5,60 |

T = mean of the spore count on non-exposed carriers (non-exposed to the Nocolyse/Nocospray System)

N = mean of the spore count on exposed carriers placed in different locations in the four rooms

R. Log = Reduction of the strain population by the disinfection system in comparison to the mean of the non exposed carriers (T)

9. CONCLUSIONS

A 6 log reduction was achieved for all 4 rooms and all 5 tested carriers for each pathogen tested.

For the ATCC 43598 strain, there was a reduction from $3 \cdot 10^7$ to $\pm 7.6 \cdot 10^1$ (76 spores remaining) which represents a log reduction of 5.78.

For the NAP1A/027 strain, there was a reduction from $1 \cdot 10^7$ to $\pm 3.9 \cdot 10^1$ (39 spores remaining) which represents a log reduction of 5.59.

For the ATCC 43598 strain, 99.99974% of the spores were eliminated (76 spores remaining out of the 30 000 000 spores). For the NAP1A/027 strain, 99.99961% of the spores were eliminated (39 spores remaining out of the 10 000 000 spores).

In order to accommodate the study, different volumes than the ones stated in the standard were used. The final number of cells was equal to the one stated in the ASTM standard.


 J.M. HUBRECHTS, Ph.D,
 Médecin hygiéniste

HIS-IZZ - Site J. BRACOPS
 Dr J.M. HUBRECHTS
 Chef du laboratoire HIS
 1 03234 71 860

J. VAN BROËK ✓
 Responsable de l'étude



10. ANNEXE 1 : DESCRIPTION OF THE TEST ROOMS (310-311-304-305)

The four test rooms were identical (same dimensions and layout) and all located on the same floor (Bracops Hospital, B Building, 3rd floor):

Standard hospital room with 1 bed, night table, chair and table, small fridge and build in cupboard.

Height of the room: 2.70 m Volume of the room: 48,8 m³



The temperature and humidity of the room were recorded during the trial as well as the distances between the appliances in the room and the carriers.

Washroom :



Seperate washroom : 6,6 m³




Hallway of the room :




Interior passageway: During the trial, the door towards the washroom was left open.

11. ANNEXE 2 : PLAN OF THE TEST ROOMS (310-311-304-305)



Location of the Carriers in the Room



Legend :

- (1) left corner of the room beside the bed
- (2) right corner of the room beside the bed
- (3) in the middle of the room, on the bed
- (4) beginning of the corridor in the room
- (5) end of the corridor in the room
- (6) Nocospray position (at a height of 40 cm)

12. ANNEXE 3: TESTING PARAMETERS

| PARAMETERS | ROOM 304 | ROOM 305 | ROOM 310 | ROOM 311 |
|--|--------------------|--------------------|--------------------|--------------------|
| At time 0 (Before the treatment) : | | | | |
| Time when the parameters were taken | 12 :30 | 12 :40 | 11 :45 | 11 :55 |
| Temperature | 23.7°C | 24.1°C | 22°C | 23.5°C |
| Humidity | 46% | 50% | 49% | 53% |
| Volume of the Room | 50 m ³ | 50 m ³ | 50 m ³ | 50 m ³ |
| Tested Product | Nocolyse | Nocolyse | Nocolyse | Nocolyse |
| Lot Number | 060509 OS | 220709 OS | 060509 OS | 220709 OS |
| Expiry Date | 05/2011 | 07/2011 | 05/2011 | 07/2011 |
| Time at the beginning of the treatment | 12 :30 | 12 :40 | 11 :45 | 11 :55 |
| Time at the end of the treatment | 12 :52 | 13 :02 | 12 :08 | 12 :17 |
| Programmed volume | 350 m ³ | 350 m ³ | 350 m ³ | 350 m ³ |
| Volume of the product used | 350 ml | 350 ml | 350 ml | 350 ml |
| At time 1 (After treatment with the Nocolyse/Nocospray System) : | | | | |
| Temperature | 25.2°C | 25.4°C | 25.6°C | 26.5°C |
| Humidity | 47% | 46% | 48% | 47% |

References

ASTM 2197:2002. Standard quantitative disk carrier test method for determining the bactericidal, virucidal, fungicidal, mycobactericidal and sporocidal activities of liquid chemical germicides. Doc