



ASAHI Frontline Guide Wires with SION TECC



	TIP LOAD (g)	TIP FLEXIBILITY	TIP LUBRICITY	SUPPORT	TORQUE RESPONSE
ASAHI SION blue	0.5	[Bar]	[Bar]	[Bar]	[Bar]
ASAHI SION	0.7	[Bar]	[Bar]	[Bar]	[Bar]
ASAHI SION black	0.8	[Bar]	[Bar]	[Bar]	[Bar]
Fielder XT-R	0.6	[Bar]	[Bar]	[Bar]	[Bar]




ASAHI SION blue




CORE	WIRE TYPE	COATINGS	WIRE O.D.	TIP LOAD
Composite SION TECC	Coil	Hydrophilic Hydrophobic Tip	0.014" (0.36mm)	0.5g

ASAHI SION




CORE	WIRE TYPE	COATING	WIRE O.D.	TIP LOAD
Composite SION TECC	Coil	Hydrophilic	0.014" (0.36mm)	0.7g

ASAHI SION black



CORE	WIRE TYPE / COATING	WIRE O.D.	TIP LOAD
Composite SION TECC	Polymer Jacket / Hydrophilic	0.014" (0.36mm)	0.8g

Fielder XT-R



CORE	WIRE TYPE / COATING	WIRE O.D.	TIP LOAD
Composite SION TECC	Polymer Jacket / Hydrophilic	0.014" (0.36mm) 0.010" (0.26mm) Tip	0.6g

• Hydrophilic = Coated with SLIP-COAT® coating

"ASAHI SION" and "Fielder" are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries.