

ASAHI

Peripheral Guide Wire

New Guide Wires

ASAHI **Gladius** ASAHI **Halberd**
ASAHI **Gaia PV** **Astato XS 40**


ASAHI's unique guide wire technology with Composite Core

SION TECC 
Leading to the NEXT


ASAHI

Peripheral Guide Wire

Exceptional guide wires for precise interventions.



**Regalia
Treasure
Astato**



ASAHI's unique guide wire technology with Composite Core
SION TECC 
Leading to the NEXT

ASAHI **Gladius**
ASAHI **Halberd**
ASAHI **Gaia PV**

Frontline

Regalia XS 1.0

0.36mm (0.014inch) Sub-total/Tortuous



Support ● ● ○ ○ ○
Torque ● ○ ○ ○ ○
Penetration ● ○ ○ ○ ○


Balanced lubricity and trackability for BTK procedures.

Tip Load
1.0 G

Polymer Jacket

Treasure Floppy

0.46mm (0.018inch) Workhorse



Support ● ● ○ ○ ○
Torque ● ● ○ ○ ○
Penetration ● ● ○ ○ ○

Precise handling with a flexible tip for tortuous anatomy.

Tip Load
4.0 G

Hybrid Coating
Silicone/Hydrophilic

Specialty

ASAHI Halberd

0.36mm (0.014inch)

Complex Lesion



Support ● ● ● ○ ○
Torque ● ● ● ● ○
Penetration ● ● ● ● ●

Composite core improves torque for navigation, and the micro-cone tip improves penetration.

Tip Load
12.0 G

Hydrophilic Coating
Uncoated Ball Tip

WITH ACT ONE

Micro-cone Tip

Treasure 12

0.46mm (0.018inch)

Complex Lesion



Support ● ● ○ ○ ○
Torque ● ● ● ● ○
Penetration ● ● ○ ○ ○

Consistent response and tactile feedback for complex lesion navigation.

Tip Load
12.0 G

Hydrophilic Coating
Uncoated Ball Tip

Astato XS 20

0.36mm (0.014inch)

Highly Calcified Lesion



Support ● ● ○ ○ ○
Torque ● ● ● ● ○
Penetration ● ● ● ● ●

Tapered tip provides excellent penetration.

Tip Load
20.0 G

Hydrophilic Coating
Uncoated Ball Tip

Tapered

Astato 30

0.46mm (0.018inch)

Highly Calcified Lesion



Support ● ● ● ○ ○
Torque ● ● ● ● ●
Penetration ● ● ● ● ●

Tapered tip provides excellent penetration.

Tip Load
30.0 G

Hydrophilic Coating
Uncoated Ball Tip

Tapered

ASAHI Gladius 0.36mm (0.014inch)

Workhorse



Support	● ● ○ ○ ○	Composite core improves tip durability and torque to treat lesions in multiple vessels.
Torque	● ● ● ○ ○	
Penetration	● ● ○ ○ ○	

Tip Load
3.0 G

Polymer Jacket



ASAHI Gladius 0.46mm (0.018inch)

Workhorse



Support	● ● ● ● ●	Composite core improves torque response with increased support for device delivery.
Torque	● ● ● ○ ○	
Penetration	● ● ○ ○ ○	

Tip Load
4.0 G

Polymer Jacket



ASAHI Halberd 0.46mm (0.018inch)

Complex Lesion



Support	● ● ● ● ○	Composite core provides precise torque for navigation with more support and micro-cone tip for proximal cap penetration.
Torque	● ● ● ● ●	
Penetration	● ● ● ● ○	


Tip Load
12.0 G

Hydrophilic Coating
Uncoated Ball Tip



ASAHI Gaia PV 0.46mm (0.018inch)

Complex Lesion



Support	● ● ● ● ○	Composite core with a lower tip load improves navigation and torque. The tapered coil and micro-cone tip makes it easier to advance.
Torque	● ● ● ● ●	
Penetration	● ● ● ● ○	

Tip Load
7.5 G

Hydrophilic Coating
Uncoated Ball Tip



Astato XS 40 0.36mm (0.014inch)

Highly Calcified Lesion




Support	● ● ○ ○ ○	Tapered tip provides excellent penetration.
Torque	● ● ● ● ○	
Penetration	★ ★ ★ ★ ★	

Tip Load
40.0 G

Hydrophilic Coating
Uncoated Ball Tip



Wire Selection Chart



Ordering Information

Product	Diameter	Catalog No.	Length	SION TECC	Tip Shape
ASAHI Gladius	0.36mm (0.014inch)	PPW14R100S	200cm	SION TECC	Straight
		PPW14R200S	235cm		
		PPW14R300S	300cm		
	0.36mm (0.014inch)	PPW14R100P	200cm		Pre-shaped
		PPW14R200P	235cm		
		PPW14R300P	300cm		
	0.46mm (0.018inch)	PPW18R100S	200cm		Straight
		PPW18R200S	235cm		
		PPW18R300S	300cm		
ASAHI Halberd	0.36mm (0.014inch)	PHW14R101S	200cm	SION TECC	Straight
		PHW14R201S	235cm		
		PHW14R301S	300cm		
	0.36mm (0.014inch)	PHW14R101P	200cm		Pre-shaped
		PHW14R201P	235cm		
		PHW14R301P	300cm		
	0.46mm (0.018inch)	PHW18R101S	200cm		Straight
		PHW18R201S	235cm		
		PHW18R301S	300cm		
	0.46mm (0.018inch)	PHW18R101P	200cm		Pre-shaped
		PHW18R201P	235cm		
		PHW18R301P	300cm		
ASAHI Gaia PV	0.46mm (0.018inch)	PHW18R102S	200cm	SION TECC	Straight
		PHW18R202S	235cm		
		PHW18R302S	300cm		
	0.46mm (0.018inch)	PHW18R102P	200cm		Pre-shaped
		PHW18R202P	235cm		
		PHW18R302P	300cm		
Regalia XS 1.0	0.36mm (0.014inch)	PAGP140000	180cm		Straight
		PAGP140300	300cm		
Astato XS 20	0.36mm (0.014inch)	PAGH143092	180cm		
		PAGH143392	300cm		
Astato XS 40	0.36mm (0.014inch)	PAGHW143094	200cm		
		PAGHW143394	300cm		
Astato 30	0.46mm (0.018inch)	PAGH18M071	180cm		
		PAGH18M371	300cm		
Treasure Floppy	0.46mm (0.018inch)	PAGH18M072	190cm		
		PAGH18M372	300cm		
Treasure 12	0.46mm (0.018inch)	PAGH18M070	180cm		
		PAGH18M370	300cm		

Indications for Use

This Product is intended to facilitate the placement and exchange of diagnostic and therapeutic devices during intravascular procedures. This device is intended for peripheral vascular use only.

Not all shapes and sizes are available in all markets. Please check availability with your sales representative. Coated with SLIP-COAT® coating.

CE 0344

© 2016 ASAHI INTECC CO., LTD.

"ASAHI", "Astato", "ASAHI Gaia", "ASAHI Gladius" and "Treasure" are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries.

Manufactured by:

ASAHI INTECC

ASAHI INTECC CO.,LTD.

3-100 Akatsuki-cho, Seto, Aichi 489-0071 Japan

Contact phone number: +81-561-86-9101

<http://www.asahi-intecc.com>

Distributed by:

Brage
MEDICAL