


ASAHI Frontline Wires



ASAHI SION blue
Well Balanced Workhorse




TIP LOAD

0.5gf



ASAHI SION
Tracks Tortuous Vessels



0.7gf




ASAHI SION black
Sliding Performance



0.8gf




Fielder XT-R
Navigates Small Vessels
and Tight Lesions



0.6gf



Fielder XT-A
Penetrates Tapered
Proximal Caps



1.0gf




ULTIMATEbros 3
Controlled Drilling



3.0gf



ASAHI Gaia First
ASAHI Gaia Second
ASAHI Gaia Third
Torque and Tip Flexibility



1.7gf
3.5gf
4.5gf



CONFIANZA PRO
Conquest Pro
CONFIANZA PRO 12
Conquest Pro 12
Penetration Technique




9.0gf
12gf

ASAHI CTO Wires


Support



ASAHI SION blue ES
Support



0.5gf



Grand Slam
Support



0.7gf

Coatings Key:  Hydrophilic Coating  Hydrophobic (Silicone) Coating  Polymer Jacket & Hydrophilic Coating

Frontline Wires



CTO Wires

