

# NOCOSPRAY

The light and portable aerial surface bio-disinfection device.



■ ■ Manufactured in France by Oxy'Pharm according to ISO 13485.

**Nocospray** is a light and portable aerial surface bio-disinfection device. It ensures aerial surface disinfection or disinsectisation when used along with our Nocolyse or Oxyphy ranges.

A dry fog allows a uniform treatment of every cm<sup>2</sup> of surfaces, as a preventive or curative measure.

Treatment of volumes from 10 to 1,000 m<sup>3</sup>.



### EASY TO USE

Easy adjustment for the volume of the room to be treated thanks to its touchpad.



### ECONOMICAL

Very small amount of product required.



### GUARANTEED TIME SAVING

No preparation required before treatment (except for cleaning) or after treatment (no need to wipe or ventilate).

## COMPATIBILITY

- Disinfectant bottles (1 L): Nocolyse, Nocolyse+, Nocolyse Food, Nocolyse One Shot.
- Insecticide bottles (1 L): Oxyphy, Oxyphy+.

## PRECAUTIONS FOR USE

Read all product information before use. Instructions for use are available on request by email: [info@oxypharm.net](mailto:info@oxypharm.net)

## TECHNICAL CHARACTERISTICS

Dimensions (L x W x H cm)	49.4 x 30 x 34.2 cm
Empty weight	5.9 kg
Packaging	
- Dimensions (L x W x H cm)	60 x 42 x 42 cm
- Weight	9 kg
Maximal treatment volume	1,000 m <sup>3</sup>
Turbine speed	22,000 rpm
Output speed of the mist	80 m/s
Average liquid flow rate	1,000 ml/h
Storage conditions	
- Temperature	0°C to 70°C
- Hygrometry	10% to 95%
Power	1,000 W
Average particle size	<5 μ

## REFERENCES

DEVICE	REFERENCE
Nocospray 230 V	2000.001
Nocospray 115 V	2000.011