




NeuroSlider®
Microcatheter DLC



Superior torqueability and pushability
Smooth and safe device delivery
Long-lasting tip shape retention


ORDERING INFORMATION | NeuroSlider® DLC

Product Name	Reference Number	ID (Inch)	OD dist. / prox. (French)	Usable Length (cm)	Tip Shape	Tip Marker
NeuroSlider® 17 DLC	01-000282	0.0165	1.9 / 2.3	155	Straight (shapeable)	2
	01-000283	0.0165	1.9 / 2.3	160	Straight (shapeable)	2
	01-000284	0.0165	1.9 / 2.3	167	Straight (shapeable)	2
NeuroSlider® 21 DLC	01-000292	0.021	2.2 / 2.6	155	Straight (shapeable)	2
	01-000293	0.021	2.2 / 2.6	160	Straight (shapeable)	2
	01-000294	0.021	2.2 / 2.6	167	Straight (shapeable)	2
NeuroSlider® 27 DLC	01-000276	0.027	3.0 / 3.1	155	Straight (shapeable)	1
NeuroSlider® 27 DLC pro*	01-000277	0.027	3.0 / 3.1	155	Straight (shapeable)	1
NeuroSlider® 39 DLC	01-000262	0.039	4.0 / 4.1	125	Straight (shapeable)	1
	01-000263	0.039	4.0 / 4.1	135	Straight (shapeable)	1
	01-000264	0.039	4.0 / 4.1	145	Straight (shapeable)	1

* For availability and device compatibility please contact your local representative from Acandis®

www.3we.de 2107

NeuroSlider® Microcatheter DLC



All changes or modifications, may they be technical or other, or changes in the availability of products are expressively reserved.

CE 0297