

ORDERING INFORMATION | ACCERO®

Labelled ACCERO® Dimensions (mm)	Reference Number	Stent Diameter (mm)	Stent Length (mm)	Recommended Vessel Diameter (mm)	Required / Recommended Microcatheters for Delivery (Inch)
2.5 × 10	01-000800	2.5	10	1.5 – 2.5	0.0165 – 0.017 NeuroSlider® 17 DLC
2.5 × 15	01-000801	2.5	15	1.5 – 2.5	
2.5 × 20	01-000802	2.5	20	1.5 – 2.5	
3.0 × 10	01-000803	3.0	10	2.0 – 3.0	
3.0 × 15	01-000804	3.0	15	2.0 – 3.0	
3.0 × 20	01-000805	3.0	20	2.0 – 3.0	
3.5 × 10	01-000806	3.5	10	2.5 – 3.5	
3.5 × 15	01-000807	3.5	15	2.5 – 3.5	
3.5 × 20	01-000808	3.5	20	2.5 – 3.5	
3.5 × 25	01-000841	3.5	25	2.5 – 3.5	
4.0 × 15	01-000810	4.0	15	3.0 – 4.0	
4.0 × 20	01-000811	4.0	20	3.0 – 4.0	
4.0 × 25	01-000845	4.0	25	3.0 – 4.0	