

Ordering information

Product	Catalog No.	O.D.	Coating	Usable length	Coil length	Radiopaque length	Tip shape
ASAHI CHIKAI 008	WAIN-CKI-008-200	0.20 / 0.26 mm (0.008 / 0.010 inch)	Hydrophilic coating (SLIP-COAT®) 180 cm	200 cm	9 cm	9 cm	Straight
ASAHI CHIKAI 10	WAIN-CKI-10-200	0.26 mm (0.010 inch)	Hydrophilic coating (SLIP-COAT®) 170 cm	200 cm	9.5 cm	3 cm	Straight
	WAIN-CKI-10-300	0.26 mm (0.010 inch)	Hydrophilic coating (SLIP-COAT®) 170 cm	300 cm	9.5 cm	3 cm	Straight
ASAHI CHIKAI	WAIN-CKI-200	0.36 mm (0.014 inch)	Hydrophilic coating (SLIP-COAT®) 170 cm	200 cm	30 cm	5 cm	Straight
	WAIN-CKI-200-RC	0.36 mm (0.014 inch)	Hydrophilic coating (SLIP-COAT®) 170 cm	200 cm	30 cm	5 cm	 Round curve
	WAIN-CKI-300	0.36 mm (0.014 inch)	Hydrophilic coating (SLIP-COAT®) 170 cm	300 cm	30 cm	5 cm	Straight
	WAIN-CKI-300-RC	0.36 mm (0.014 inch)	Hydrophilic coating (SLIP-COAT®) 170 cm	300 cm	30 cm	5 cm	 Round curve
ASAHI CHIKAI black	WAIN-CKI-200-BS	0.36 mm (0.014 inch)	Polymer jacket 150 cm / Hydrophilic coating (SLIP-COAT®) 170 cm	200 cm	30 cm	3 cm	 Round curve
	WAIN-CKI-200-BA	0.36 mm (0.014 inch)	Polymer jacket 150 cm / Hydrophilic coating (SLIP-COAT®) 170 cm	200 cm	30 cm	3 cm	 Angled 90°
ASAHI CHIKAI black 14 soft tip	AIN-CKI-200-B-SFT	0.36 mm (0.014 inch)	Polymer jacket 150 cm / Hydrophilic coating (SLIP-COAT®) 170 cm	200 cm	30 cm	3 cm	 Round curve
ASAHI CHIKAI black 18	WAIN-CKI-18-200-BS	0.45 / 0.36 mm (0.018 / 0.014 inch)	Polymer jacket 150 cm / Hydrophilic coating (SLIP-COAT®) 170 cm	200 cm	34 cm	5 cm	 Round curve
ASAHI CHIKAI black 18 soft tip	AIN-CKI-18-200-SFT	0.45 / 0.36 mm (0.018 / 0.014 inch)	Polymer jacket 150 cm / Hydrophilic coating (SLIP-COAT®) 170 cm	200 cm	34 cm	5 cm	 Round curve

Product	Catalog No.	O.D.	Coating	Usable Length	Coil Length	Radiopaque Length	Tip Shape
ASAHI CHIKAI X 010	AIN-CKX-10-200-R	0.26mm (0.010inch)	Hydrophilic coating (SLIP-COAT®) 180cm	200cm	16cm	16cm	Re-shapeable round curve