

ASAHI INTECC


ASAHI Frontline and CTO Guide Wires

Category	Product Name	Description	Tip Load
ASAHI Frontline Wires	ASAHI SION blue	Well Balanced Workhorse	0.5g
	ASAHI SION	Tracks Tortuous Vessels	0.7g
	NEW ASAHI SION black	Sliding Performance	0.8g
	NEW Fielder XT-R	Navigates Small Vessels and Tight Lesions	0.6g
ASAHI CTO Wires	NEW Fielder XT-A	Penetrates Tapered Proximal Caps	1.0g
	ULTIMATEbros 3	Controlled Drilling	3.0g
	NEW ASAHI Gaia First ASAHI Gaia Second ASAHI Gaia Third	Torque and Tip Flexibility	1.7g, 3.5g, 4.5g
	CONFIANZA PRO Conquest Pro	Penetration Technique	9.0g
	CONFIANZA PRO 12 Conquest Pro 12		12.0g
	Support	Grand Slam Support	

ASAHI Frontline Wires



ASAHI CTO Wires



Coatings Key

- Hydrophilic Coating
- Hydrophobic (Silicone) Coating
- Polymer Jacket & Hydrophilic Coating