



Ordering Information

ASAHI Masters PARKWAY SOFT




Product	Catalog No.	Usable Length	Length	O.D. of Distal shaft / O.D. of Proximal shaft	Shaft I.D.	Max. Guide Wire	Hydrophilic Coating Length	Tip Shape	Max. Pressure
ASAHI Masters PARKWAY SOFT	WMST105-18PWSF	105cm	112cm	0.66mm/1.98Fr 0.94mm/2.8Fr	0.55mm/0.022inch	0.46mm/0.018inch	65cm	Straight	6895kPa/1000psi
	WMST125-18PWSF	125cm	132cm	0.66mm/1.98Fr 0.94mm/2.8Fr	0.55mm/0.022inch	0.46mm/0.018inch	65cm	Straight	6895kPa/1000psi
	WMST150-18PWSF	150cm	157cm	0.66mm/1.98Fr 0.94mm/2.8Fr	0.55mm/0.022inch	0.46mm/0.018inch	100cm	Straight	6895kPa/1000psi
	WMST45A-18PWSSF	105cm	112cm	0.66mm/1.98Fr 0.94mm/2.8Fr	0.55mm/0.022inch	0.46mm/0.018inch	65cm	Angled	6895kPa/1000psi
	WMST45A-18PWSF	125cm	132cm	0.66mm/1.98Fr 0.94mm/2.8Fr	0.55mm/0.022inch	0.46mm/0.018inch	65cm	Angled	6895kPa/1000psi

ASAHI Masters PARKWAY HF



Product	Catalog No.	Usable Length	Length	O.D. of Distal shaft / O.D. of Proximal shaft	Shaft I.D.	Max. Guide Wire	Hydrophilic Coating Length	Tip Shape	Max. Pressure
ASAHI Masters PARKWAY HF	WMST105-27HF	105cm	112cm	0.88mm/2.6Fr 0.93mm/2.8Fr	0.69mm/0.027inch	0.53mm/0.021inch	65cm	Straight	6895kPa/1000psi
	WMST125-27HF	125cm	132cm	0.88mm/2.6Fr 0.93mm/2.8Fr	0.69mm/0.027inch	0.53mm/0.021inch	65cm	Straight	6895kPa/1000psi
	WMST150-27HF	150cm	157cm	0.88mm/2.6Fr 0.93mm/2.8Fr	0.69mm/0.027inch	0.53mm/0.021inch	100cm	Straight	6895kPa/1000psi

ASAHI Masters PARKWAY HF KIT



Product	Catalog No.	Micro Catheter Usable Length	Micro Catheter Length	O.D. of Distal shaft / O.D. of Proximal shaft	Shaft I.D.	Guide Wire Shaft O.D.	Guide Wire Length	Tip shape of GuideWire	Max. Pressure
ASAHI Masters PARKWAY HF KIT	WMST105-27HFK	105cm	112cm	0.88mm/2.6Fr 0.93mm/2.8Fr	0.69mm/0.027inch	0.53mm/0.021inch	140cm	Angled 80°	6895kPa/1000psi
	WMST125-27HFK	125cm	132cm	0.88mm/2.6Fr 0.93mm/2.8Fr	0.69mm/0.027inch	0.53mm/0.021inch	160cm	Angled 80°	6895kPa/1000psi
	WMST105-27HFKM	105cm	112cm	0.88mm/2.6Fr 0.93mm/2.8Fr	0.69mm/0.027inch	0.53mm/0.021inch	140cm	Multi Curve	6895kPa/1000psi
	WMST125-27HFKM	125cm	132cm	0.88mm/2.6Fr 0.93mm/2.8Fr	0.69mm/0.027inch	0.53mm/0.021inch	160cm	Multi Curve	6895kPa/1000psi



ASAHI and MastersPARKWAY are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries.

Manufactured by:



ASAHI INTECC CO.,LTD.
http://www.asahi-intecc.com

Distributed by:

