

# ASPIRATION CATHETER

For thrombus-containing lesions & SVGs

## STENTYS AC




- ~26% increase in extraction speed<sup>1</sup>
- Largest effective lumen area for extraction<sup>2</sup>
- High kink-resistance



STENTYS

# STENTYS® Aspiration Catheter

Ordering Code	Guiding Catheter Compatibility	Product Specification
AC01-6-01	6F	Hydrophilic Stiff



- A - Atraumatic soft distal tip
- B - Guidewire compatibility: 0.014" (0.36mm)
- C - Outer diameter: 0.067" (1.70mm)
- D - Proximal shaft for excellent pushability

STENTYS® AC is indicated for the containment and aspiration of embolic material (thrombus/debris) while performing PTCA and/or stenting procedures, including in saphenous vein grafts.

Extraction lumen area : 0.95mm<sup>2</sup>  
Aspiration speed : 1.6 cc/s

The STENTYS® AC system is CE marked.

1. Compared with Medtronic Export XT
2. Compared with Terumo Eliminate™, Medtronic Export XT & AP, Medtronic (Invatec) Diver CE

CAUTION: Indications, precautions, warnings and instructions for use can be found in the product labelling supplied with the device.  
Not available for sale in the United States.  
© STENTYS S.A. 2012

ST-AC-120521

STENTYS S.A.  
25, rue de Choiseul  
75002 Paris  
France  
Tel: +33 1 44 53 99 42  
Fax: +33 1 44 53 99 24

STENTYS Inc.  
Suite #109  
103 Carnegie Center  
Princeton, NJ 08540  
Tel: +1 (609) 853 0100  
Fax: +1 (609) 275 6155

[www.stentys.com](http://www.stentys.com)