

# ASAHI SION blue ES

P T C A G U I D E W I R E

# Supportive Front

## ASAHI SION blue ES PTCA GUIDE WIRE


An extra support type wire joins the line-up of ASAHI SION Series.  
A wire that stabilizes the guide catheter for stress-free device delivery.

### Characteristics

- **Contribute to a stable device delivery**  
Well balanced reinforced shaft enhances guiding catheter disengagement during device delivery.
- **Inherits the maneuverability of SION TECC**  
Usable as a frontline guide wire with easy control for branch selection and minimal loss of steerability in tortuous vessels.

### Support

Stable device delivery by improving shaft reinforce.



All the data were obtained by company standardized test, which may differ from industry standardized tests.  
All the data do not guarantee that all devices have exactly the same performance with the samples used for tests.



# line guidewire

## SION TECC

- Enhanced tip durability
- Combination of flexibility and torque force
- Reduced whipping of the wire due to superior torque response

ASAHI's unique guide wire technology with Composite Core



“Small Changes, Great Possibilities” **SION TECC**   
Leading to the NEXT

Variety of choices that suits you best

**ASAHI SION series guidewire**



# Structure & Ordering Information



Employs the ASAHI brand rope coil, which provides torqueability, flexibility and resiliency.

Product	Catalog No.	Diameter	Usable Length	Coil Length	Radiopaque Length	Tip Shape	Tip Load
ASAHI SION blue ES	AHW14R017S	0.36mm (0.014inch)	190cm	8.5cm	3cm	Straight	0.5gf
	AHW14R017J	0.36mm (0.014inch)	190cm	8.5cm	3cm	J	0.5gf
	AHW14R317S	0.36mm (0.014inch)	300cm	8.5cm	3cm	Straight	0.5gf
	AHW14R317J	0.36mm (0.014inch)	300cm	8.5cm	3cm	J	0.5gf

CE 0344

Not all shapes and sizes are available in all markets. Please check availability with your sales representative.  
Coated with SLIP-COAT® coating.

© 2016 ASAHI INTECC CO., LTD.  
\*ASAHI\*, \*SION\* and \*ACT ONE\* are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries.

Manufactured by:

Distributed by:

**ASAHI INTECC**

**ASAHI INTECC CO.,LTD.**

3-100 Akatsuki-cho, Seto, Aichi 489-0071 Japan

Contact phone number: +81-561-86-9101

<http://www.asahi-intecc.com>