

Ordering Information

Product	Catalog No. 6.5Fr	Catalog No. 7.5Fr	Catalog No. 8.5Fr	Shape Code	Usable Length (cm)
Judkins Left	G6JL35-0-L100 G6JL40-0-L100 G6JL45-0-L100 G6JL50-0-L100	G7JL35-0-L100 G7JL40-0-L100 G7JL45-0-L100 G7JL50-0-L100	G8JL35-0-L100 G8JL40-0-L100	JL3.5 JL4.0 JL4.5 JL5.0	100 100 100 100
Judkins Left (Short Tip Type)	G6JL35-0-L100S G6JL40-0-L100S G6JL45-0-L100S G6JL50-0-L100S	G7JL35-0-L100S G7JL40-0-L100S G7JL45-0-L100S G7JL50-0-L100S	G8JL35-0-L100S G8JL40-0-L100S	JL3.5 ST JL4.0 ST JL4.5 ST JL5.0 ST	100 100 100 100
Judkins Right	G6JR35-0-L100 G6JR40-0-L100 G6JR45-0-L100 G6JR50-0-L100	G7JR35-0-L100 G7JR40-0-L100 G7JR45-0-L100 G7JR50-0-L100	G8JR35-0-L100 G8JR40-0-L100	JR3.5 JR4.0 JR4.5 JR5.0	100 100 100 100
Judkins Right (Short Tip Type)	G6JR35-0-L100S G6JR40-0-L100S G6JR45-0-L100S G6JR50-0-L100S	G7JR35-0-L100S G7JR40-0-L100S G7JR45-0-L100S G7JR50-0-L100S	G8JR35-0-L100S G8JR40-0-L100S	JR3.5 ST JR4.0 ST JR4.5 ST JR5.0 ST	100 100 100 100
Amplatz Left	G6AL07-0-L100 G6AL10-0-L100 G6AL15-0-L100 G6AL20-0-L100	G7AL07-0-L100 G7AL10-0-L100 G7AL15-0-L100 G7AL20-0-L100	G8AL10-0-L100	AL0.75 AL1.0 AL1.5 AL2.0	100 100 100 100
Amplatz Left (Short Tip Type)	G6AL07-0-L100S G6AL10-0-L100S G6AL15-0-L100S G6AL20-0-L100S	G7AL07-0-L100S G7AL10-0-L100S G7AL15-0-L100S G7AL20-0-L100S		AL0.75 ST AL1.0 ST AL1.5 ST AL2.0 ST	100 100 100 100
Special Curve	G6SC35-0-L100 G6SC40-0-L100	G7SC35-0-L100 G7SC40-0-L100		SC3.5 SC4.0	100 100
Hockey Stick	G6HS00-0-L100	G7HS00-0-L100		HS	100
Power Backup	G6PB30-0-L100 G6PB35-0-L100 G6PB40-0-L100	G7PB30-0-L100 G7PB35-0-L100 G7PB40-0-L100		PB3.0 PB3.5 PB4.0	100 100 100
Super Power Backup	G6SP30-0-L100 G6SP35-0-L100 G6SP40-0-L100	G7SP30-0-L100 G7SP35-0-L100 G7SP40-0-L100	G8SP35-0-L100 G8SP40-0-L100	SP3.0 SP3.5 SP4.0	100 100 100
Multipurpose	G6MP10-0-L100 G6MP20-0-L100	G7MP10-0-L100 G7MP20-0-L100	G8MP10-0-L100 G8MP20-0-L100	MP1.0 MP2.0	100 100

Not all shapes and sizes are available in all markets. Please check availability with your sales representative.

Manufactured by:



ASAHI INTECC CO.,LTD.

3-100 Akatsuki-cho, Seto, Aichi 489-0071 Japan
TEL: +81-561-48-5551 FAX: +81-561-48-5552
<http://www.asahi-intecc.com>

Distributed by:



ASAHI INTECC CO., LTD. EUROPE Strawinskylaan 967 (WTC Tower D-9th floor) 1077 XX Amsterdam, The Netherlands TEL: +31-(0)20-794-0640 FAX: +31-(0)20-794-0641	ASAHI INTECC USA, INC. 2500 Red Hill Avenue, Suite 210 Santa Ana, California, 92705, USA TEL: +1-949-756-8901 FAX: +1-949-756-8165	ASAHI INTECC (THAILAND) CO., LTD. 158/1 Moo5, Bangkadi Industrial Park, Tiwanon Road, Tambol Bangkadi, Amphur Muang, Pathumthani 12000, Thailand TEL: +66-2-501-1302 FAX: +66-2-501-1305	ASAHI INTECC CO., LTD. SINGAPORE BRANCH 3 Bishan Place #06-01 Bishan CPE Building Singapore 579838 TEL: +65-6250-4248 FAX: +65-6250-4247	ASAHI INTECC CO., LTD. BEIJING REPRESENTATIVE OFFICE Room 18, 10th Floor, Twin Towers (East), B12 Jianguomenwai Ave., Beijing, China 100022 TEL: +86-10-5123-5100 FAX: +86-10-5123-5177	ASAHI INTECC CO., LTD. HONG KONG BRANCH Unit 501-2, 5/F, Hon Kwok Jordan Centre No 7 Hillwood Road, Tsim Sha Tsui, Kowloon Hong Kong TEL: +852-2316-2510 FAX: +852-2316-2485
--	--	---	--	--	---

PTCA GUIDING CATHETER

SheathLess
Eaucath



The Minimally Invasive TRI Approach – Without Sacrificing the Space You Need

PTCA GUIDING CATHETER


SheathLess
Eaucath

Sheathless Eaucath represents a new concept in guide catheters. Perform a full range of procedures – from the simplest to the most complex – via the radial approach, without sacrificing function or lumen size.

Features

- Substantial Outer Layer & Large Inner Lumen**
Enhanced kink resistance and backup support with a large inner lumen for easy device delivery.
- Hydrophilic Coating**
Enhanced catheter trackability even in tortuous vessels, and reduced incidence of spasm.
- ASAHI Braiding**
Two different braiding patterns provide optimal torque and flexibility, while helping to maintain the lumen.

6.5Fr : 0.070"
7.5Fr : 0.081"
8.5Fr : 0.090"



Smaller punctures for larger catheters – good for the patient, good for the physician.

The OD of a Sheathless Eaucath catheter is similar to the OD of a sheath introducer 2 sizes smaller. Especially useful for ad-hoc PCI situations.

