

E-mail mottatt. 23.05.2019 – kl 11:46

Dear Sven

Virkon contains 10-15% of the surfactant, & Perasafe contains < 3%.

We consider that Perasafe is however the better cleaning agent on hard surfaces overall, but neither of the products are promoted as cleaners (they are however both considered & sold as disinfectants - Virkon as a biocide, under the Biocidal Products Regulation (BPR, Regulation (EU) 528/2012) , & Perasafe as a Medical Device under ISO 13485:2016, ISO 9001:2015).

Best regards

Mark Squire
Technical Marketing Manager Disinfection
Antec International Limited
Windham Road
Sudbury
CO10 2XD
United Kingdom

Antec - a company of the LanXess Group

Email mark.squire@lanxess.com
Phone +44 1787 377305