

SIMPLY UNIQUE

ICAD* Treatment with
NeuroSpeed® PTA Balloon Catheter and CREDO® Stent



- Approved for acute (LVO) and elective stenting
- Repositionable up to 90% deployment
- New diameter (5.0 mm) and longer lengths (25 mm, 30 mm)

* ICAD Intracranial Atherosclerotic Disease

xcandis®

ENGINEERING STROKE SOLUTIONS

ORDERING INFORMATION

Labelled CREDO® Dimensions (mm)	Reference Number	Stent Diameter (mm)	Stent Length (mm)	Recommended Vessel Diameter (mm)	Required Catheter for Delivery
3.0 × 15	01-000930	3.0	15	2.0–2.5	NeuroSpeed® PTA Balloon Catheter
3.0 × 20	01-000931	3.0	20	2.0–2.5	NeuroSpeed® PTA Balloon Catheter
3.0 × 25	01-000932	3.0	25	2.0 - 2.5	NeuroSpeed® PTA Balloon Catheter
3.0 × 30	01-000933	3.0	30	2.0 - 2.5	NeuroSpeed® PTA Balloon Catheter
4.0 × 15	01-000940	4.0	15	2.5–3.5	NeuroSpeed® PTA Balloon Catheter
4.0 × 20	01-000941	4.0	20	2.5–3.5	NeuroSpeed® PTA Balloon Catheter
4.0 × 25	01-000942	4.0	25	2.5 - 3.5	NeuroSpeed® PTA Balloon Catheter
4.0 × 30	01-000943	4.0	30	2.5 - 3.5	NeuroSpeed® PTA Balloon Catheter
5.0 × 15	01-000950	5.0	15	3.5–4.5	NeuroSpeed® PTA Balloon Catheter
5.0 × 20	01-000951	5.0	20	3.5–4.5	NeuroSpeed® PTA Balloon Catheter
5.0 × 25	01-000952	5.0	25	3.5–4.5	NeuroSpeed® PTA Balloon Catheter
5.0 × 30	01-000953	5.0	30	3.5–4.5	NeuroSpeed® PTA Balloon Catheter

All sizes feature HRF (High Radial Force)

Labelled NeuroSpeed® Dimensions (mm)	Reference Number	Balloon Diameter (mm)	Balloon Working Length (mm)	ID (Inch)	OD dist. / prox. (French)	Usable Length (cm)
1.5 × 8	01-000605	1.5	8	0.0165	2.7 / 3.7	150
2.0 × 8	01-000600	2.0	8	0.0165	2.7 / 3.7	150
2.5 × 8	01-000601	2.5	8	0.0165	2.7 / 3.7	150
3.0 × 8	01-000602	3.0	8	0.0165	2.7 / 3.7	150
3.5 × 8	01-000603	3.5	8	0.0165	2.7 / 3.7	150
4.0 × 8	01-000604	4.0	8	0.0165	2.7 / 3.7	150

Inflation Pressure (bar)	NeuroSpeed® Diameter (mm)					
	1.5	2.0	2.5	3.0	3.5	4.0
2.0	1.21	1.72	2.09	2.42	3.06	3.26
4.0	1.37	1.84	2.33	2.78	3.25	3.72
6.0	1.50*	2.00*	2.50*	3.00*	3.50*	4.00*
8.0	1.67	2.16	2.65	3.22	3.69	4.23
10.0	1.85	2.27	2.75	3.38	3.83	4.37
12.0	2.02	2.39	2.87	3.54	3.97**	4.53**
14.0	2.20**	2.52**	2.98**	3.73**	–	–

* Nominal pressure ** Rated burst pressure

All changes or modifications, may they be technical or other, or changes in the availability of products are expressly reserved.

Distributed by:



ACANDIS GmbH
Theodor-Fahrner-Str. 6
75177 Pforzheim
Germany

Tel: +49 7231 155 00 0
Fax: +49 7231 155 00 129
E-Mail: info@acandis.com
www.acandis.com