

Varenavn:

ASAHI wires

Produktgruppe:

Perifere Guide Wires

Produsent:

ASAHI



Varenavn	Varenummer (code)	Produktbeskrivelse	Wire lengde (cm)	HydrofilCoating g	Polymer sleeve	Spring Coil lengde (cm)	Tupp form	Tupp "vekt"	Röntgentett lengde (cm)
ASAHI Gladius 0.014"	PPW14R100S	ASAHI Gladius 0.014" x 200 cm rett tupp	200	50 cm	40 cm	11	rett	3 g	3
	PPW14R200S	ASAHI Gladius 0.014" x 235 cm rett tupp	235	50 cm	40 cm	11	rett	3 g	3
	PPW14R300S	ASAHI Gladius 0.014" x 300 cm rett tupp	300	50 cm	40 cm	11	rett	3 g	3
	PPW14R100P	ASAHI Gladius 0.014" x 200 cm J-bøy	200	50 cm	40 cm	11	J	3 g	3
	PPW14R200P	ASAHI Gladius 0.014" x 235 cm J-bøy	235	50 cm	40 cm	11	J	3 g	3
	PPW14R300P	ASAHI Gladius 0.014" x 300 cm J-bøy	300	50 cm	40 cm	11	J	3 g	3
0.018"	PPW18R100S	ASAHI Gladius 0.018" x 200 cm rett tupp	200	10 cm	10 cm	5	rett	4 g	4
	PPW18R200S	ASAHI Gladius 0.018" x 235 cm rett tupp	235	10 cm	10 cm	5	rett	4 g	4
	PPW18R300S	ASAHI Gladius 0.018" x 300 cm rett tupp	300	10 cm	10 cm	5	rett	4 g	4
Gladus MG 14 PV 0.014"	PP14R003P	Asahi Gladius MG 14 PV 0.014" x 190 cm J-bøy	190	41 cm	41 cm	8,5	J	3 g	3
	PP14R203P	Asahi Gladius MG 14 PV 0.014" x 235 cm J-bøy	235	41 cm	41 cm	8,5	J	3 g	3
	PP143303P	Asahi Gladius MG 14 PV 0.014" x 300 cm J-bøy	300	41 cm	41 cm	8,5	J	3 g	3
Gladus MG 14 PV ES 0.014"	PP14R004P	Asahi Gladius MG 14 PV ES 0.014" x 190 cm J-bøy	190	10 cm	10 cm	3	J	3 g	3
	PP14R204P	Asahi Gladius MG 14 PV ES 0.014" x 235 cm J-bøy	235	10 cm	10 cm	3	J	3 g	3
	PP143304P	Asahi Gladius MG 14 PV ES 0.014" x 300 cm J-bøy	300	10 cm	10 cm	3	J	3 g	3
Gladus MG 18 PV ES 0.018"	PP18R004P	Asahi Gladius MG 18 PV ES 0.014" x 190 cm J-bøy	190	10 cm	10 cm	4,5	J	3 g	4,5
	PP18R204P	Asahi Gladius MG 18 PV ES 0.014" x 235 cm J-bøy	235	10 cm	10 cm	4,5	J	3 g	4,5
	PP18R304P	Asahi Gladius MG 18 PV ES 0.014" x 300 cm J-bøy	300	10 cm	10 cm	4,5	J	3 g	4,5
ASAHI Gaia PV 0.018"	PHW18R102S	ASAHI Gaia PV 0.018" x 200 cm rett tupp	200	16 cm	-	15	rett	7,5 g	15
	PHW18R202S	ASAHI Gaia PV 0.018" x 235 cm rett tupp	235	16 cm	-	15	rett	7,5 g	15
	PHW18R302S	ASAHI Gaia PV 0.018" x 300 cm rett tupp	300	16 cm	-	15	rett	7,5 g	15
	PHW18R102P	ASAHI Gaia PV 0.018" x 200 cm J-bøy	200	16 cm	-	15	J	7,5 g	15
	PHW18R202P	ASAHI Gaia PV 0.018" x 200 cm J-bøy	235	16 cm	-	15	J	7,5 g	15
	PHW18R302P	ASAHI Gaia PV 0.018" x 300 cm J-bøy	300	16 cm	-	15	J	7,5 g	15
ASAHI Halberd 0.014"	PHW14R101S	ASAHI Halberd 0.014" x 200 cm rett tupp	200	33 cm	-	5	rett	12 g	5
	PHW14R201S	ASAHI Halberd 0.014" x 235 cm rett tupp	235	33 cm	-	5	rett	12 g	5
	PHW14R301S	ASAHI Halberd 0.014" x 300 cm rett tupp	300	33 cm	-	5	rett	12 g	5
	PHW14R101P	ASAHI Halberd 0.014" x 200 cm J-bøy	200	33 cm	-	5	J	12 g	5
	PHW14R201P	ASAHI Halberd 0.014" x 235 cm J-bøy	235	33 cm	-	5	J	12 g	5
	PHW14R301P	ASAHI Halberd 0.014" x 300 cm J-bøy	300	33 cm	-	5	J	12 g	5
0.018"	PHW18R101S	ASAHI Halberd 0.018" x 200 cm rett tupp	200	12 cm	-	11	rett	12 g	11
	PHW18R201S	ASAHI Halberd 0.018" x 235 cm rett tupp	235	12 cm	-	11	rett	12 g	11
	PHW18R301S	ASAHI Halberd 0.018" x 300 cm rett tupp	300	12 cm	-	11	rett	12 g	11
	PHW18R101P	ASAHI Halberd 0.018" x 200 cm J-bøy	200	12 cm	-	11	J	12 g	11
	PHW18R201P	ASAHI Halberd 0.018" x 235 cm J-bøy	235	12 cm	-	11	J	12 g	11
	PHW18R301P	ASAHI Halberd 0.018" x 300 cm J-bøy	300	12 cm	-	11	J	12 g	11



Act One Technology: enhances tip durability and reduces whipping of the wire due to superior torque control.

Varenavn:

ASAHI wires

Produktgruppe:

Perifere Guide Wires,

Produsent:

ASAHI



Varenavn	Varenummer (code)	Produktbeskrivelse	Wire lengde (cm)	HydrofilCoatin g	Polymer sleeve	Spring Coil lengde (cm)	tupp form	tupp "vekt"	Røntgentett lengde (cm)
Regalia XS 1.0 0.014"	PAGP140000	Regalia XS 1,0 - 0.014" x 180 cm rett tupp	180	22 cm	22 cm	12	rett	2 g	3
	PAGP140300	Regalia XS 1,0 - 0.014" x 300 cm rett tupp	300	22 cm	22 cm	12	rett	2 g	3
Astato XS 20 0.014"	PAGH143092	Astato XS 20 - 0.014" x 180 cm rett tupp	180	17 cm	-	17	rett	20 g	17
	PAGH143392	Astato XS 20 - 0.014" x 300 cm rett tupp	300	17 cm	-	17	rett	20 g	17
Astato XS 40 0.014"	PAGHW143094	Astato XS 40 - 0.014" x 200 cm rett tupp	200	17 cm	-	17	rett	40 g	17
	PAGHW143394	Astato XS 40 - 0.014" x 300 cm rett tupp	300	17 cm	-	17	rett	40 g	17
Astato 30 0.018"	PAGH18M071	Astato 30 - 0.018" x 180 cm rett tupp	180	15 cm	-	15	rett	30 g	15
	PAGH18M371	Astato 30 - 0.018" x 300 cm rett tupp	300	15 cm	-	15	rett	30 g	15
Treasure Floppy 0.018"	PAGH18M072	Treasure Floppy 0.018" x 190 cm rett tipp	190	15 cm	-	15	rett	4 g	3
	PAGH18M372	Treasure Floppy 0.018" x 300 cm rett tipp	300	15 cm	-	15	rett	4 g	3
Treasure 12 0.018"	PAGH18M070	Treasure 12 - 0.018" x 180 cm rett tipp	180	15 cm	-	15	rett	12 g	15
	PAGH18M370	Treasure 12 - 0.018" x 300 cm rett tipp	300	15 cm	-	15	rett	12 g	15
Sion black 0.014"	APW14R010S	Asahi SION black 190 cm x 0.014 rett tupp	190	40 cm	20 cm	12	rett	0,8 g	3
	APW14R310S	Asahi SION black 300 cm x 0.014 rett tupp	300	40 cm	20 cm	12	rett	0,8 g	3
	APW14R010J	Asahi SION black 190 cm x 0.014 J-bøy	190	40 cm	20 cm	12	J	0,8 g	3
	APW14R310J	Asahi SION black 300 cm x 0.014 J-bøy	300	40 cm	20 cm	12	J	0,8 g	3



Act One Technology: enhances tip durability and reduces whipping of the wire due to superior torque control.