

Vascular access system

y Connector




Unique push-pull system
Clean, safe and fast



Committed to excellence in developing Arterial Therapeutic Systems

INTENDED USE

To maintain haemostasis before and during the introduction, the use and withdrawal of diagnostic/interventional devices outer diameter $\leq 7F$ (2.33mm / 0.092") during interventional procedures.

**REFERENCES**

Reference	Internal reference	Description	Insertion tool ⁽¹⁾	Torquer ⁽²⁾	Plastic peelable insertion tool	Packed/Bulk	Sterile/Non-sterile
18020005	KESU0001	Y connector 20 PSI	x	-	-	P	S
18020006	KESU0003	Y connector 50 PSI	x	-	-	P	S
18020007	KESU0004	Y connector 00 PSI	x	-	x	P	S
18020008	KESU0007	Y connector 20 PS	-	-	-	P	S
18020009	KESU0008	Y connector 50 PS	-	-	-	P	S
18020010	KESU0009	Y connector 00 PS	-	-	x	P	S
18020011	KESU0002	Y connector 20 PSTI	x	x	-	P	S
18020012	KESU0005	Y connector 50 PSTI	x	x	-	P	S
18020013	KESU0006	Y connector 00 PSTI	x	x	-	P	S

⁽¹⁾ Compatible with guidewires up to 0.018"


⁽²⁾ Compatible with guidewires 0.010" to 0.035"

Important information: Prior to use, refer to the "instructions for use" supplied with these devices for indications, contraindications, side effects, suggested procedures, warnings and precautions. For healthcare professionals only. Class IIa medical device.

CONTACT US >>>

 qco.net/qmd/arthesys/

 info@arthesys.com

 +33(0)141118777

 4 rue René Razel
91400 Saclay
France