

The fast pressure bandage for radial puncture site care with/without transparent inspection possibility



Indication

Premofix[®] uno Radialis for compression of radial puncture sites after invasive diagnosis and therapy.

Characteristics

- ready to use bandage for left or right
- integrated, transparent compression segment (Radialis T)
- with special transparent thumb-pressure-plate (Radialis T)
- visual control already during application
- easy and quick to apply
- no slipping of the backstrap through the eye
- stable and secure by unique material combination
- tear-resistant even during strong extension
- bacteria-impermeable material but kind to skin
- high transverse stability
- especially soft strap material
- optimizes and accelerates the treatment of the patient
- enables an earlier mobilization of the patient

Premofix[®] uno Radialis – for left or right

Product	Description	Art.No.	Shipping unit
Premofix [®] uno Radialis C	Single-use – pressure bandage to apply pressure on punctured radial arteries	DV1079C	40 pieces

