




ACCLINO® flex plus
Stent



High visibility
More treatment options
Compatible with 0.0165" – 0.027" ID microcatheters

ORDERING INFORMATION | ACCLINO® flex plus

Labelled ACCLINO® flex plus Dimensions (mm)	Reference Number	Stent Diameter (mm)	Stent Length (mm)	Recommended Vessel Diameter (mm)	Required / Recommended Microcatheters for Delivery (Inch)
[HRF] 3.0 × 15 3.0 × 20 3.0 × 25 3.0 × 30 3.0 × 35	01-001122	3.0	15	1.5 – 2.5	0.0165 – 0.017 NeuroSlider® 17 NeuroSlider® 17 DLC
	01-001123	3.0	20	1.5 – 2.5	
	01-001124	3.0	25	1.5 – 2.5	
	01-001125	3.0	30	1.5 – 2.5	
	01-001126	3.0	35	1.5 – 2.5	
[HRF] 3.5 × 15 3.5 × 20 3.5 × 25 3.5 × 30 3.5 × 35	01-001132	3.5	15	1.5 – 3.0	
	01-001133	3.5	20	1.5 – 3.0	
	01-001134	3.5	25	1.5 – 3.0	
	01-001135	3.5	30	1.5 – 3.0	
	01-001136	3.5	35	1.5 – 3.0	
[HRF] 4.0 × 15 4.0 × 20 4.0 × 25 4.0 × 30 4.0 × 35	01-001142	4.0	15	2.5 – 3.5	
	01-001143	4.0	20	2.5 – 3.5	
	01-001144	4.0	25	2.5 – 3.5	
	01-001145	4.0	30	2.5 – 3.5	
	01-001146	4.0	35	2.5 – 3.5	
[HRF] 4.5 × 15 4.5 × 20 4.5 × 25 4.5 × 30 4.5 × 35	01-001152	4.5	15	2.5 – 4.0	
	01-001153	4.5	20	2.5 – 4.0	
	01-001154	4.5	25	2.5 – 4.0	
	01-001155	4.5	30	2.5 – 4.0	
	01-001156	4.5	35	2.5 – 4.0	
[HRF] 5.0 × 15 5.0 × 20 5.0 × 25 5.0 × 30 5.0 × 35	01-001162	5.0	15	3.0 – 4.5	
	01-001163	5.0	20	3.0 – 4.5	
	01-001164	5.0	25	3.0 – 4.5	
	01-001165	5.0	30	3.0 – 4.5	
	01-001166	5.0	35	3.0 – 4.5	
[HRF] 5.5 × 20 5.5 × 25 5.5 × 30 5.5 × 35	01-001173	5.5	20	3.5 – 5.0	
	01-001174	5.5	25	3.5 – 5.0	
	01-001175	5.5	30	3.5 – 5.0	
	01-001176	5.5	35	3.5 – 5.0	
[HRF] 6.5 × 20 6.5 × 25 6.5 × 30 6.5 × 35	01-001193	6.5	20	4.0 – 6.0	0.021 NeuroSlider® 21 NeuroSlider® 21 DLC
	01-001194	6.5	25	4.0 – 6.0	
	01-001195	6.5	30	4.0 – 6.0	
	01-001196	6.5	35	4.0 – 6.0	
[HRF] 8.0 × 20* 8.0 × 30* 8.0 × 40* 8.0 × 60*	01-001213	8.0	20	6.0 – 7.0	0.027 NeuroSlider® 27 (DLC)
	01-001215	8.0	30	6.0 – 7.0	
	01-001217	8.0	40	6.0 – 7.0	
	01-001221	8.0	60	6.0 – 7.0	

HRF: High Radial Force – compared to ACCLINO® flex plus Stents within the same recommended vessel diameter.

* For availability please contact your local representative from Acandis®

All changes or modifications, may they be technical or other, or changes in the availability of products are expressly reserved.

CE 0297