



NeuroBridge[®]
Catheter




Precise navigation
Excellent stability and support
Soft rounded pre-shaped tip (25° MP)

ORDERING INFORMATION | NeuroBridge®

Product Name	Reference Number	ID (Inch)	OD dist. (French / Inch)	OD prox. (French / Inch)	Usable Length (cm)	Total Length (cm)	Tip Shape
NeuroBridge® 39	01-000508	0.039	3.9/0.051	4.2/0.055	125	131	Multi-Purpose 25°
	01-000509	0.039	3.9/0.051	4.2/0.055	135	141	Multi-Purpose 25°
	01-000510	0.039	3.9/0.051	4.2/0.055	145	151	Multi-Purpose 25°
NeuroBridge® 52	01-000511	0.052	5.0/0.066	5.3/0.070	115	121	Multi-Purpose 25°
	01-000512	0.052	5.0/0.066	5.3/0.070	125	131	Multi-Purpose 25°
	01-000513	0.052	5.0/0.066	5.3/0.070	135	141	Multi-Purpose 25°
NeuroBridge® 65	01-000514	0.065	6.1/0.080	6.3/0.083	115	121	Multi-Purpose 25°
	01-000515	0.065	6.1/0.080	6.3/0.083	125	131	Multi-Purpose 25°

www.3we.de 2107



All changes or modifications, may they be technical or other, or changes in the availability of products are expressly reserved.

CE 0297

09-000061 EN 06/2021